

=> d his

(FILE 'HOME' ENTERED AT 10:52:00 ON 28 JUN 2001)
FILE 'CA' ENTERED AT 10:52:06 ON 28 JUN 2001
L1 4258 S CLHB OR BLOOD(2A)SUBSTITUT? OR HBOC OR(CROSSLINK? OR CROSS LINK?)
(2A)(HB OR HEMOGLOB? OR HAEMOGLOB?)
L2 77 S L1 AND INTERFER?
L3 39 S L2 NOT PY>1996
L4 5 S L2 NOT L3 AND PATENT/DT
L5 44 S L3-4
FILE 'MEDLINE' ENTERED AT 10:55:43 ON 28 JUN 2001
L6 26 S L3
FILE 'BIOSIS' ENTERED AT 10:56:00 ON 28 JUN 2001
L7 19 S L3
FILE 'CA, MEDLINE, BIOSIS' ENTERED AT 10:56:28 ON 28 JUN 2001
L8 66 DUP REM L5 L6 L7 (23 DUPLICATES REMOVED)

=> d 18 bib,ab 1-66

L8 ANSWER 9 OF 66 MEDLINE
AN 96409702 MEDLINE
TI Selective determination of chloride and bromide ions in serum by cyclic voltammetry.
AU Arai K; Kusu F; Noguchi N; Takamura K; Osawa H
CS School of Pharmacy, Tokyo University of Pharmacy and Life Science, 1432-1 Horinouchi, Hachioji, Tokyo, 192-03, Japan.
SO ANALYTICAL BIOCHEMISTRY, (1996 Aug 15) 240 (1) 109-13.
AB Voltammetric determination of chloride and bromide ions in serum was made using a one-body type Ag electrode, whose surface was covered with a dialysis membrane for preventing interference from proteins. This determination was based on measurement of the charge of the reduction wave of silver halide formed on the Ag electrode surface in a halide ion solution during a cathodic potential sweep. In this method, depending on the potential sweep range, the charge of the reduction wave of AgI, AgI + AgBr, or AgI + AgBr + AgCl can be measured. Here, the method was used to determine the chloride and bromide ions in artificial serum. Linear concentration ranges were 1.0×10^{-2} - 2.4×10^{-1} M for chloride ions and 5.0×10^{-3} - 1.8×10^{-2} M for bromide ions. The correlation coefficient was 0.999 in each case. The relative standard deviation for chloride ion (1.0×10^{-1} M) was 1.9% (n = 10), and that for bromide ion (1.5×10^{-2} M) was 2.5% (n = 10). Bromide ions at a relatively high concentration as the component of serum had no effect on the determination of the chloride ion in the serum. This method is thus useful for the selective determination of chloride and bromide ions in serum.

L8 ANSWER 11 OF 66 MEDLINE
AN 96101281 MEDLINE
TI Artificial icteric plasmas: unreliable indicators for interference with creatinine assay on Beckman CX3.
AU van't Sant P; Kreutzer H J
CS Lab. of Clin. Chem., Bosch Medicentrum, 's-Hertogenbosch, The Netherlands.
SO CLINICAL CHEMISTRY, (1995 Dec) 41 (12 Pt 1) 1773-4.

L8 ANSWER 19 OF 66 BIOSIS COPYRIGHT 2001 BIOSIS
AN 1993:403890 BIOSIS
TI Modification of clinical chemistry methods to overcome interferences from diaspirin crosslinked hemoglobin (DCLHb.
AU Leissing, Nancy; Mattia-Goldberg, Cynthia; Oskroba, Diana

CS Baxter Healthcare Corp., Round Lake, IL 60073 USA
SO Clinical Chemistry, (1993) Vol. 39, No. 6, pp. 1144. Meeting Info.: 45th
National Meeting of the American Association for Clinical Chemistry, Inc.
New York, New York, USA July 11-15, 1993 ISSN: 0009-9147.

L8 ANSWER 20 OF 66 BIOSIS COPYRIGHT 2001 BIOSIS
AN 1993:403884 BIOSIS
TI Modification of hematology methods to overcome interferences from diaspirin
crosslinked hemoglobin (DCLHb).
AU Leissing, Nancy; Mattia-Goldberg, Cynthia; Balma, Diane
CS Healthcare Corp., Round Lake, IL 60073 USA
SO Clinical Chemistry, (1993) Vol. 39, No. 6, pp. 1142. Meeting Info.: 45th
National Meeting of the American Association for Clinical Chemistry, Inc.
New York, New York, USA July 11-15, 1993 ISSN: 0009-9147.

L8 ANSWER 21 OF 66 BIOSIS COPYRIGHT 2001 BIOSIS
AN 1993:403822 BIOSIS
TI Hemoglobin based blood substitutes: Interference with routine chemical
tests.
AU Christenson, Robert H.; Gregory, Linda C.; Duh, Show-Hong; MacKenzie,
Colin; Trump, Benjamin F.
CS Dep. Pathol. Anesth., University Maryland Sch. Med., Baltimore, MD 21201
USA
SO Clinical Chemistry, (1993) Vol. 39, No. 6, pp. 1129. Meeting Info.: 45th
National Meeting of the American Association for Clinical Chemistry, Inc.
New York, New York, USA July 11-15, 1993 ISSN: 0009-9147.

L8 ANSWER 32 OF 66 CA COPYRIGHT 2001 ACS
AN 113:146775 CA
TI Interference in endotoxin and fungal polysaccharide assays from blood
products and antimicrobial agents
AU Ikemura, Katsumi; Tanaka, Hiroshi; Yoshioka, Toshiharu; Sugimoto, Tsuyoshi
CS Med. Sch., Osaka Univ., Osaka, 553, Japan
SO Rinsho Byori (1990), 38(1), 87-92
AB A conventional chromogenic Limulus test (CCLT, normal plasma value: <10
pg/mL) detects both of endotoxin and β -1,3-glucan, while a chromogenic
endotoxin-specific test (EST: <3 pg/mL) detects only endotoxin. Two
sources of pseudo-pos. result in CCLT and/or EST were studied; (1) factor
G-activating material found in blood products, which is indicative of
fungal infections, and (2) interaction in diazo-coupling process of the
tests by antibiotics or sulfonamides. CCLT and EST values obtained from 13
vials of blood products including human albumin, antithrombin III, and
antihemophilic globulin were 284.6 pg/mL and 20.2 pg/mL, resp. High CCLT
values accompanied by low EST values indicate factor G-activating material
in the blood products rather than endotoxin. Sample blanks of 32 plasma
specimens obtained from 19 patients were 23.7 pg/mL by CCLT and 19.2 pg/mL
by EST. All the patients were administered sulfamethoxazole/trimethoprim
(SMX/TMP), aztreonam (AZT), ceftazidime (CAZ), imipenem/cilastatin
(IPM/CS), or silver sulfadiazine. SMX/TMP, AZT, CAZ, IPM/CS, and
cefmenoxime were demonstrated to react with the diazo-coupling reagents
used in CCLT and EST, and to cause false-pos. result in the tests. Thus,
the investigation of blood product and sample blank is required for the
assay of endotoxins or fungal polysaccharides using CCLT and EST.

L8 ANSWER 49 OF 66 MEDLINE
AN 84020957 MEDLINE
TI Effects of Fluosol-DA (artificial blood) on clinical chemistry tests and
instruments.

RB1, 193
AU Mullins R E; Hutton P S; Conn R B
SO AMERICAN JOURNAL OF CLINICAL PATHOLOGY, (1983 Oct) 80 (4) 478-83.
AB Artificial blood must be added to the list of therapeutic agents that produce interference with diagnostic laboratory tests. Fluosol-DA (Alpha Therapeutic Corp., Los Angeles, CA), a stable 20% emulsion of perfluorocarbons in aqueous medium, is being evaluated in clinical trials as a blood substitute in the United States. We investigated its effects in blood and serum samples on test results and instruments in the clinical chemistry laboratory. The 20% emulsion was added to blood or serum specimens in amounts corresponding to the replacement of in-vivo plasma volumes of 10-50%, concentrations that would be expected in blood samples obtained from patients who have received Fluosol. Observed interferences mimicked those caused by high triglyceride concentrations in serum specimens: interference with chemical reactions and generation of spurious absorbance readings because of turbidity. These types of errors are often additive, and the cumulative effect may cause either erroneously high or low values for the analytes concerned. Because Fluosol may be used widely, although infrequently, for patients refusing blood transfusions on religious grounds and for patients with rare antibodies to red blood cells who require transfusion, laboratories analyzing specimens containing Fluosol should be aware of the potential errors.

L8 ANSWER 62 OF 66 CA COPYRIGHT 2001 ACS
AN 70:84799 CA
TI Effects of drugs on laboratory values
AU Lubran, Myer
CS Univ. of Chicago, Chicago, Ill., USA
SO Med. Clin. N. Amer. (1969), 53(1), 211-22
AB This review consists primarily of listings of drugs affecting clin. lab. results, together with minimal comments relative to the means whereby these effects are achieved. Included are drugs interfering with test procedures, with the activity of anticoagulants, with thyroid function tests, with liver function tests, drugs inducing liver microsomal enzymes, drugs affecting porphyrin metabolism, drugs inducing blood dyscrasias, drugs acting on the central nervous system, analgesics, autonomic drugs, hypotensive agents, diuretics, chemotherapeutic and antibiotic agents, antifungal and antiprotozoal drugs, blood and blood substitutes, gastrointestinal drugs, vitamins, hormones and substitutes, cholesterol-lowering drugs, and misc. 5 references.

=> log y

STN INTERNATIONAL LOGOFF AT 10:57:50 ON 28 JUN 2001